

# ULTIMATE YEAST DENATURING IP

---

Where possible, we prefer to use denaturing IPs to recover labeled proteins from pulse-chase experiments. Use of denaturing IPs not only reduces background, but also makes protein harvest easy—yeast are simply boiled in an SDS-containing solution to liberate total proteins. This is also much safer than bead beating radioactive yeast samples.

---

## A. EXTRACTING YEAST PROTEINS.

Prepare: 250  $\mu$ l SHD per set of four samples (ie., one pulse-chase).

1. Resuspend each frozen yeast pellet in 50  $\mu$ l SHD by pipeting. Heat at 95°C for 10 minutes.
2. Spin briefly. Transfer supernatant to a fresh snap-lock tube.

## B. IMMUNOPRECIPITATION.

Prepare: (i) 1.2 ml ice-cold TNN with inhibitors per sample (5 ml for 4 samples).  
(ii) 2 sets of labeled tubes.  
(iii) 25  $\mu$ l of Protein A Sepharose combined with 25  $\mu$ l of Protein G sepharose per sample (125  $\mu$ l of each for 4 samples). Wash two times in TNN and resuspend in 25  $\mu$ l of TNN (125  $\mu$ l for 4 samples).

1. Add 200  $\mu$ l of TNN to each sample, vortex, and spin for 2 minutes in the cold.
2. Transfer supernatant to a fresh tube containing 1 ml of TNN plus 50  $\mu$ l of Protein A/G sepharose.
3. Rotate in cold room for 1 hr.
4. Spin lysates in microfuge for 20 seconds in the cold. Transfer supernatant to a fresh tube.
5. Add antibody as appropriate (we use 0.5  $\mu$ l of 12CA5 ascites for HA-tagged proteins).
6. Rotate in cold room for 1 hour.

Prepare: (i) 25  $\mu$ l of Protein A Sepharose combined with 25  $\mu$ l of Protein G sepharose per sample (125  $\mu$ l of each for 4 samples). Wash two times in TNN and resuspend in 25  $\mu$ l of TNN (125  $\mu$ l for 4 samples).

7. Spin samples briefly. Add 50  $\mu$ l of Protein A/G Sepharose mix. Rotate in cold room for 1 hr.

Prepare: (i) 3 ml TNN with Complete Inhibitor per sample (15 ml for 4 samples).  
(ii) 100  $\mu$ l of 2xLaemmli buffer with 50 mM DTT (NB: this is one-quarter the usual amount of DTT).

8. Spin samples for 2 minutes on setting number 3 in microfuge. Remove supernatant.

9. Wash sepharose beads by adding 1 ml TNN and rocking back and forth 10 times.

10. Spin samples for 2 minutes on setting number 3. Remove supernatant. Repeat wash two more times.

11. Aspirate all liquid from the beads using a 26G needle. Resuspend in 100  $\mu$ l 2 x Laemmli. Boil and load 50  $\mu$ l on SDS-PAGE.

**C. SOLUTIONS.**TNN

50 mM Tris (pH 7.5)  
250 mM NaCl  
5 mM EDTA  
0.5% NP-40

**For 500 ml**

25 ml 1M Tris (pH 7.5)  
25 ml 5 M NaCl  
5 ml 0.5 M EDTA  
25 ml 10 % NP-40

**Add fresh each time:**

0.4 mg/ml Pefabloc  
10 µg/ml Leupeptin  
10 µg/ml Pepstatin  
5 µg/ml Aprotinin

1:100 of 40 mg/ml stock  
1:1000 of 10 mg/ml stock  
1:500 of 2 mg/ml stock  
1:2000 of 10 mg/ml stock

or....Add 1 COMPLETE INHIBITOR TABLET per 25-30 ml for IP washes.

Note: We use individual protease inhibitors when performing the extraction and initial IP steps. For the washes, we use Complete protease inhibitor tablets (Roche: 1 tablet per 25 ml TNN).

SHD

2% SDS  
90 mM HEPES (pH 7.5)

**For 1 ml**

200 µl 10% SDS  
90 µl 1M HEPES (pH7.5)

**Add fresh each time:**

30 mM DTT  
0.4 mg/ml Pefabloc

30 µl 1M DTT  
1:100 of 40 mg/ml stock