

ULTIMATE HIS-UB ASSAY FOR YEAST

This assay is performed to detect ubiquitylated proteins in yeast. Yeast are first transformed with the vector pUB221 (from Dan Finley: Daniel_Finley@hms.harvard.edu), which expresses polyhistidine-tagged Ub under the control of the copper-inducible CUP1 promoter. Yeast are grown, induced with copper, and harvested. Total ubiquitylated proteins are then recovered by nickel-affinity chromatography, and specific proteins probed for by Western blotting.

A. GROWTH AND INDUCTION.

1. Grow 200 ml of yeast culture overnight. When OD_{600} is 0.4–1.0, induce culture with 400 μ l 250 mM $CuSO_4$. Grow 4 hr at 30°C.
2. Measure OD_{600} . Collect 5 ml for total lysate and 195 ml for purification of His-tagged proteins. Spin down yeast, remove all liquid, and freeze pellets in liquid nitrogen. Store pellets at $-80^\circ C$ until ready.

B. PURIFICATION.

His-tagged proteins.

1. To each large (195 ml pellet) frozen sample, add Buffer A to give a final OD of approximately 150.
2. Vortex. Add glass beads so there is approximately 1-2 mm of free liquid above the beads.
3. Vortex each sample: 20s vortex/30 s ice. Repeat 6 times.
4. Add 2 ml Buffer A. Vortex.
5. Spin 3000 rpm bench top centrifuge for 15 minutes at 4°C.
6. Transfer supernatant to fresh tube.
7. Measure the [protein] in 1 μ l of remaining sample, using the Bradford Assay.

8. For each sample, transfer an equal amount of protein to a fresh tube. Add 250 μ l of 50 % Ni-NTA beads (Quiagen; equilibrated in Buffer A).
9. Rotate on bench for 1-2 hr.
10. Pellet beads by spinning 2 minutes at 3000 rpm in microfuge.
11. Aspirate supernatant. Wash beads in 1 ml Buffer A.
12. Wash beads for a total of three times in Buffer A, three times in Buffer A/TI, and one time in buffer TI.
13. Pellet beads after final wash. Remove all liquid from beads.
14. Resuspend in 200 μ l of 2 x Laemmli with 0.2 M imidazole. Boil for five minutes. Run 100 μ l on appropriate SDS-PAGE gel.

Total Proteins.

1. Bring tubes up from -80 degree freezer on ice.
2. To each tube, add appropriate amount of EZ Buffer™. Mix by pipeting up and down.
3. Place in 95 degree block for 10 minutes.
4. Transfer 1 μ l of each sample to a fresh tube. Measure protein concentration using Bradford assay. Include 1 μ l of EZ Buffer™ in each of the standards. Freeze the remaining samples until ready to run gel.
5. Transfer equal amounts of protein for each sample (~75 μ g) to a fresh tube. Add 2 x Laemmli buffer to make the final volume approximately 60 μ l. Run on SDS-PAGE gel.

For both the His-tagged and the total proteins, probe for the presence of your proteins by Western blotting.

C. SOLUTIONS.

Buffer A

6 M guanidine-HCl
0.1 M Na₂HPO₄/NaH₂PO₄
10 mM imidazole
Adjust pH to 8.0 with NaOH.

For 250 ml

143.3 g guanidine HCl
25 ml 1M Na₂HPO₄/NaH₂PO₄.H₂O
2.5 ml 1M imidazole

TI

25 mM Tris.Cl
20 mM imidazole
Adjust pH to 6.8.

For 200 ml

5 ml 1M tris.Cl (pH 6.8)
4 ml 1M imidazole

Buffer A/TI

1 volume Buffer A
3 volumes TI

EZ Buffer™

0.06M Tris-HCl, pH 6.8
10% (v/v) glycerol
2% (w/v) SDS
5% (v/v) 2-mercaptoethanol