

ULTIMATE MAMMALIAN DENATURING IP

Where possible, we prefer to use denaturing IPs to recover labeled proteins from pulse-chase experiments. Use of denaturing IPs reduces background considerably. If, however, you prefer to use non-denaturing conditions, see the 'Ultimate Mammalian Non-Denaturing IP' protocol. The protocol here is described for four sets of pulse-chase experiments (16 samples).

A. EXTRACTING AND DENATURING PROTEINS.

Prepare: (i) 2.0 ml TSD
(ii) Labeled snap-cap microfuge tubes

1. Do not thaw cells before making extract. To frozen cell pellet, add 100 μ l of TSD. Mix up and down by pipetting. Boil 10 minutes.
2. Spin lysates in microfuge for 5 minutes at room temperature. Transfer supernatant to a fresh tube.

B. IMMUNOPRECIPITATION.

Prepare: (i) 25 ml ice-cold TNN.
(ii) Labeled snap-cap microfuge tubes
(iii) 450 μ l of Protein A Sepharose combined with 450 μ l of Protein G sepharose.
Washed 2 times in TNN and resuspended in 450 μ l TNN.

1. To each sample, add 1.2 ml of TNN and mix by pipetting. Add 50 μ l of Protein A/G sepharose mix. Rotate in the cold for 1 hr.
2. Spin lysates in microfuge at top speed for 20 seconds. Transfer supernatant to a fresh tube.
3. Add antibody (*eg.*, 0.5 μ l 12CA5 ascites)
4. Rotate in cold room for 1 hour.

Prepare: (i) 450 μ l of Protein A Sepharose combined with 450 μ l of Protein G sepharose.
Washed 2 times in TNN and resuspended in 450 μ l TNN.

5. Add 50 μ l of Protein A/G Sepharose mix. Rotate in cold room for 1 hr.

Prepare: (i) 50 ml ice-cold TNN (with Complete Inhibitor)
(ii) 3 ml 2x Laemmli buffer (with 50 mM DTT)

6. Spin samples 10 seconds in the cold (or 2 minutes on setting number 3).
Remove supernatant.

7. Wash sepharose beads by adding 1 ml TNN and rocking back and forth 10 times.

8. Spin samples 10 seconds in the cold (or 2 minutes on setting number 3).
Remove supernatant. Repeat wash two more times.

9. Aspirate all liquid from the beads using a 26G needle. Resuspend in 100 μ l of 2 x Laemmli (50 mM DTT). Boil and load 50 μ l on SDS-PAGE.

Note that the usual concentration in 2 x Laemmli is 200 mM. For these experiments, the concentration is dropped to 50 mM.

C. SOLUTIONS.TNN

50 mM Tris (pH 7.5)
250 mM NaCl
5 mM EDTA
0.5% NP-40

For 500 ml

25 ml 1M Tris (pH 7.5)
25 ml 5 M NaCl
5 ml 0.5 M EDTA
25 ml 10 % NP-40

Add fresh each time:

0.4 mg/ml Pefabloc
10 µg/ml Leupeptin
10 µg/ml Pepstatin
5 µg/ml Aprotinin

1:100 of 40 mg/ml stock
1:1000 of 10 mg/ml stock
1:500 of 2 mg/ml stock
1:2000 of 10 mg/ml stock

or....Add 1 COMPLETE INHIBITOR TABLET per 25-30 ml for IP washes.

Note: We use individual protease inhibitors when performing the extraction and initial IP steps. For the washes, we use Complete protease inhibitor tablets (Roche: 1 tablet per 25 ml TNN).

TSD

50 mM Tris (pH 7.5)
1% SDS
5 mM DTT

For 200 ml

10 ml 1M Tris (pH 7.5)
10 ml 20 % SDS
ADD 1:200 of 1M DTT FRESH