

ULTIMATE MAMMALIAN CELL PULSE-CHASE

This assay is performed to directly measure the stability of a protein in mammalian cells. The assay works well with transiently-expressed proteins, providing they have multiple-methionine residues and can be expressed at high levels in transfected cells. For proteins that can only be expressed at low levels (or proteins that have a few methionine residues) we use the methionine-fortified HAM tag (available from Bill Tansey: tansey@cshl.edu) to facilitate labeling. For detection of endogenous proteins, you may need to increase the amount of ^{35}S -labeled methionine (up to 500 μl per 6 cm dish). Labeled proteins can be recovered by denaturing immunoprecipitation (our favorite—see ‘Ultimate Mammalian Denaturing IP’ protocol) or non-denaturing immunoprecipitation (see ‘Ultimate Mammalian Non-Denaturing IP protocol’).

A. CELL SET-UP.

This protocol is designed for use with ~50% confluent 6xcm dishes of cells. If cells are to be transfected, they should be transfected 24-40 hr before the labeling begins. Labeling is done in methionine- and cysteine-free DMEM—if your cells cannot grow in DMEM, you may substitute another methionine- and cysteine-free media. Typically, each pulse-chase requires four plates of cells—one for the zero time point (pulse) and three for the subsequent time points. This protocol is described for four different pulse-chase determinations (ie., total of 16 plates).

B. PULSE-CHASE.

- Prepare:
- (i) 500 ml warm PBS.
 - (ii) 20 ml warm Met/Cys-free DMEM/5% dialyzed FCS (with Gln if needed). This DMEM and dialyzed FCS come from Gibco.
 - (iii) 900 μl of warm Tran ^{35}S -label (NEN).
 - (iv) 100 ml of regular DMEM + FCS + 2 mM Methionine + 2 mM Cysteine.
 - (v) 100 ml ice-cold PBS.
 - (vi) “Shoebox” for labeling.
 - (vii) 16 labeled microfuge tubes (on ice).
 - (viii) Disposable cell scrapers.
 - (ix) Liquid nitrogen.
 - (x) Stopwatch.
 - (xi) Disposable ^{35}S -waste (liquid and sharps).

1. Aspirate media from each plate of cells. Wash cells 2 times in 5 ml warm PBS.

2. Aspirate PBS. Add 1 ml warm DMEM (met/cys-free)/5% dialyzed FCS. Incubate cells for 30 minutes at 37°C. This is the starvation period.
3. Add 50 µl of Tran³⁵S-label (at roughly 10 mCi/ml; equals 500 µCi total) to each plate. Place cells in shoebox and incubate for 40 minutes at 37°C.
4. Carefully remove radioactive media. Wash cells on plate with 2 ml warm PBS. Remove PBS. Add 5 ml warm DMEM+FCS+2 mM Methionine+2 mM Cysteine to all plates except the zero "0" samples (5 plates total). Start timer.
5. Immediately remove all "0" plates. Return other plates to 37°C.
6. Remove media from "0" plates. Wash cells with 5 ml ice-cold PBS. Remove PBS. Scrape cells in 0.5 ml ice-cold PBS. Transfer to labeled microfuge tube and spin 45 seconds in the cold.
7. Carefully remove supernatant. Freeze cell pellet in liquid nitrogen and store at -80°C until ready to perform the immunoprecipitation.
8. Perform similar collections at different points during the chase (30, 60, 90 minutes is a good starting point).

Use the ULTIMATE MAMMALIAN CELL DENATURING IP or ULTIMATE MAMMALIAN CELL NON-DENATURING IP protocol to recover your labeled protein.

Note that the Met/Cys-free DMEM needs to have glutamine added. Also, to the chase media, add Met and Cys fresh (from 200 mM stocks in water) just before the chase.