

ULTIMATE YEAST CHIP ASSAY

This assay is for detection of proteins bound to specific regions of chromatin. There are many variations of this assay; this one works well for us.

A. GROW CELLS.

1. Inoculate 5 ml cultures of each yeast in appropriate media. Grow to saturation at appropriate temperature (usually 30°C).

2. Inoculate 50 ml cultures (250 ml flask) of selective media with 3 different amounts of the original starter culture. For cultures growing at 30°C, use 10, 30, and 100 μ l of the original starter; for cultures growing at 25°C, use 30, 100, and 300 μ l of the original culture. Grow for ~16 hr at appropriate temperature. The idea is to have at least one culture the next day that is at OD 0.5-1.0.

B. FIX CELLS.

Prepare:

25 ml of 37.5% PFA fresh by mixing:

9.4 g paraformaldehyde (stored at 4°C).

50 μ l 10N NaOH

15 ml PBS

Heat to 65°C until PFA is dissolved (7 on hotplate; takes ~20 minutes). Take to 25 ml final volume with PBS; filter through 0.45 μ M filter.

Ice-cold PBS

Benchtop centrifuge at 4 degrees

1. Select cultures that are between OD₆₀₀ 0.5 and 1.0. Place culture on rotating platform at room temperature.

2. With culture rotating, add 1.5 ml of 37% PFA dropwise to each culture. Continue to rotate for 15 minutes.

Prepare:

Labeled 50 ml falcon tubes, one for each culture

Labeled 2 ml screwcap tubes, one for each culture

3. Add 3ml of 2.5M glycine to each culture. Continue to rotate for 5 minutes at room temperature.

4. Transfer cell suspension to 50 ml conical tube. Spin top speed for five minutes at 4°C.

5. Wash cell pellet twice in 50 ml ice-cold PBS.

6. After the final wash, transfer cell pellet (in 1 ml PBS) to a 2 ml screw-cap microfuge tube. Place on ice until ready to proceed with lysis (but don't wait too long).

C. LYSE CELLS.

Prepare:

Ice-cold lysis buffer + freshly-added inhibitors (450 μ l per sample)

21G needles

Labeled siliconized microfuge tubes (no caps)

Labeled siliconized microfuge tubes (caps)

1. Pellet cells. Resuspend in 400 μ l of ice-cold LYSIS BUFFER.

2. Add 500 μ l cold, acid-washed, glass beads (500 μ M diameter). Close tube.

3. Place tubes (in pairs, so that machine is balanced) into the BioSpec bead beater, set up in the cold room. Lyse cells with 4 x 40 second pulses, with the bead beater set to 'homogenize'. Leave 30 seconds between pulses.

4. Puncture the bottom of the 2 ml screwcap tube with a hot 21G needle. Place this tube on top of a siliconized 1.5 ml eppendorf tube (with no cap), and place both into a 15 ml centrifuge tube.

5. Spin the tubes for 5 minutes at top-speed in the benchtop centrifuge.

6. Remove the top tube, and transfer the lysate (not pellet) to a new siliconized 1.5 ml eppendorf tube. Keep tube on ice.

D. SHEAR CHROMATIN.

Prepare:

Check sonicator, clean tip.

2 sets of label siliconized eppendorf tubes.

1. Keeping each tube on ice, sonicate (setting 5; microtip) for 10 seconds (count to 11), then return to ice for a minimum of 60 seconds. Repeat this procedure 2 times, for a total of 20 seconds sonication [NB: this time is calibrated for a specific sonicator, and designed to give DNA fragments

between 500 and 1 kb in length. Other crosslinking conditions and machines may require different sonication times].

2. Spin all tubes at top speed in microcentrifuge. 5 minutes at 4°C. Expect to see a small pellet. Transfer supernatants to fresh siliconized tubes. Store on ice until immunoprecipitation.

E. IMMUNOPRECIPITATION.

Prepare:

2 sets of labeled siliconized microfuge tubes; one labeled 'INPUT'

Ice-cold lysis buffer + freshly-added inhibitors (50 µl per sample)

Protein A agarose/Protein G agarose. Mix 30 µl of protein A agarose per sample with an equal volume of protein G agarose. Wash three times in lysis buffer and resuspend in 1 volume of lysis buffer (+inhibitors).

1. To each sample, add 50 µl of Protein A/G agarose. Rotate at 4°C for 1 hr.

2. Spin 2000 rpm for 5 minutes at 4°C. Transfer S/N to fresh tube.

3. Remove 50 µl from each sample to fresh tubes labeled INPUT. Freeze at -20°C until step G.

4. To each sample, add appropriate antibody. For 12CA5, use 10 μ l of ascites diluted 1:10 in lysis buffer. For 9E10, use 10 μ l of purified Ab-1 from Oncogene Sciences.

5. Rotate at 4°C for three hours.

Prepare:

Protein A agarose/Protein G agarose. Mix 30 μ l of protein A agarose per sample with an equal volume of protein G agarose. Wash three times in lysis buffer and resuspend in 1 volume of lysis buffer (+inhibitors).

6. Spin briefly. Add 50 μ l of protein A/G agarose. Incubate at 4°C for 1 hr to collect immune complexes.

F. WASHING IP.

Prepare:

*LYSIS Buffer + Complete Inhibitor. Prepare 3 ml per IP. Dissolve 1 Complete tablet per 25 ml.
Labeled siliconized microfuge tubes.*

1. Spin tubes 3000 rpm for 2 minutes. Remove supernatant. Wash beads as follows:

2 x 5 min (room temp) in 1 ml LYSIS Buffer

1 x 5 min (room temp) in 1 ml DOC Buffer

1 x 5 min (room temp) in 1 ml TE

Always use separate tip for washing!

2. At last wash, transfer pellets in TE to a fresh tube. Spin down. Remove supernatant with 26G needle. Proceed immediately to step G.

G. REVERSE CROSSLINKS.

Prepare:

Labeled siliconized microfuge tubes.

65°C waterbath

1. Add 50 μ l of TES to the beads, and incubate at 65°C for 10 minutes.
2. Spin 5 minutes at top speed in microfuge. Transfer supernatant to fresh tube.
3. Add 150 μ l of TES to the beads. Mix. Spin. Transfer supernatant to the tube in step 2.
4. Add 150 μ l of TES to the INPUT tubes.
5. Incubate all tubes at 65°C overnight.

Prepare:

2 sets of Labeled SNAPCAP microfuge tubes.

6. Spin tubes briefly. Transfer to snapcap tube containing 25 μ l of 10 mg/ml proteinase K and 200 μ l of TE. Incubate at 37°C for 2 hours.

7. Extract with 400 μ l of phenol/chloroform. Transfer supernatant to a fresh tube containing 44 μ l of 3 M NaOAc and 1 μ l of glycogen (20 mg/ml).

8. Add 1 ml ice-cold ethanol and precipitate in dry-ice ethanol bath. Pellet DNA, wash in 70 % ethanol, vacuum dry (briefly), and resuspend in 40 μ l TE containing (1:200) RNase. Store at -20°C until PCR.

SOLUTIONS

2.5 M GLYCINE

(take to pH8.0 with NaOH to get to dissolve)

ICE-COLD PBS

LYSIS BUFFER

50 mM HEPES/KOH (pH 7.5)	_____
500 mM NaCl	_____
1 mM EDTA	_____
1% Triton X-100	_____
0.1% DOC	_____
0.1% SDS	_____

PLUS INHIBITORS ADDED FRESH

0.4 mg/ml Pefabloc	1:100 of 40 mg/ml stock
10 µg/ml Leupeptin	1:1000 of 10 mg/ml stock
10 µg/ml Pepstatin	1:500 of 2 mg/ml stock
5 µg/ml Aprotinin	1:2000 of 10 mg/ml stock

TES

50 mM Tris-Cl (pH 8.0)

10 mM EDTA

1 % SDS

DOC Buffer

10 mM Tris.Cl (pH 8.0)

0.25M LiCl

0.5% NP-40

0.5% DOC

1 mM EDTA
